


Series DX1

2-sectioned recovery and patient transport table with hydraulic height adjustment

Model 1090-00/H white powder coated (RAL 9010)
Model 1090-03/H white aluminium powder coated (RAL 9006)
Model 1090-04/H grey aluminium powder coated (RAL 9007)



Both pictures:
Model 1090-04/H
with option: Storage basket
(B1850)

Basic equipment:

- Pedals for foot-controlled height positioning on both sides of the table
- Movable on twin castors, Ø 100 mm, with central lock, 1 castor with direction lock
- Push handles at head part and foot part, RAL 9006 powder coated
- 2 side guards, lateral lowerable, RAL 9006 powder coated, length: 820 mm
- 4 wall protection bumpers, Ø 200 mm
- Upholstery medium-soft (high-quality cold foam cover)
- Dewert-hygiene-equipment



Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable by gas spring (0° - 65°)
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 60 mm



Technical data:

- Manual foot operated height adjustment: 500 mm bis 950 mm by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: 700 mm, width overall: 780 mm
- Length of lying surface: 1950 mm, length overall: 2060 mm
- Load capacity: 225 kg, according to IEC 60601 tested with fourfold static load
Load capacity head part: 90 kg



Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1181
Without basic mobility



Art. no. B1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 30 mm



Art. no. B1199
Footbar for central lock, at foot end

Paper roll holder



Art. no. 011K / 013K
Paper roll holder, mounted under head part, for width 700 mm / 800 mm

Upholstery options

Art. no. 1090-016
Width 800 mm

Art. no. B1770 / B1780
Soft covering for pleasant resting comfort, for width 700 mm/ 800 mm, thickness 80 mm

Add-on elements



Art. no. B1410
1 piece of side rail, length: 350 mm



Art. no. B1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. B1450
Infusion pole, without attachment clamp



art. no B1850
Storage basket,
measures:
L. 530 x W. 370 x H. 80 mm