


**NEW**

## Series DX1


2-sectioned recovery and patient transport table  
with hydraulic height adjustment

- Model 1093-00/H white powder coated (RAL 9010)
- Model 1093-03/H white aluminium powder coated (RAL 9006)
- Model 1093-04/H grey aluminium powder coated (RAL 9007)



Both pictures:  
Model 1093-04/H  
with basic equipment;  
upholstery colour: 64-712 slate grey

### Basic equipment:

- Pedals for foot-controlled height positioning on both sides of the table
- Movable on twin castors, Ø 100 mm, with central lock, 1 castor with direction lock
- Push handles at head part and foot part, RAL 9006 powder coated
- 2 side guards, lateral lowerable, RAL 9006 powder coated, length: 820 mm
- 4 wall protection bumpers, Ø 200 mm
- Upholstery medium-soft (high-quality cold foam cover)
- Dewert-hygiene-equipment 


### Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable by gas spring (0° - +75°) (sitting position)
- Length head part: 750 mm
- Fixed foot part
- Length foot part: 1200 mm
- Upholstery thickness: 70 mm




### Technical data:

- Manual foot operated height adjustment: 540 - 990 mm by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: 700 mm, width overall: 780 mm
- Length of lying surface: 1950 mm, length overall: 2060 mm
- Load capacity: 225 kg, according to IEC 60601 tested with fourfold static load  kg
- Load capacity head part: 90 kg

### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

### Mobility



Art. no. 1181  
Without basic mobility



Art. no. B1185  
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 30 mm



Art. no. B1199  
Footbar for central lock, at foot end

### Paper roll holder



Art. no. 011K / 013K  
Paper roll holder, mounted under head part, for width 700 mm / 800 mm

### Upholstery options

Art. no. 1093-016  
Width 800 mm

Art. no. B1770 / B1780  
Soft covering for pleasant resting comfort, for width 700 mm/ 800 mm, thickness 80 mm

### Add-on elements



Art. no. B1410  
1 piece of side rail, length: 350 mm



Art. no. B1430  
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. B1450  
Infusion pole, without attachment clamp



art. no B1850  
Storage basket,  
measures:  
L. 530 x W. 370 x H. 80 mm