

Tables for echocardiography

2-sectioned table for echocardiography with fixed height

... for examination, treatment ...

Model 117-00 and 177-00 white powder coated (RAL 9010)

Model 117-01 and 177-01 chromium-plated

Model 117-03 and 177-03 white aluminium powder coated (RAL 9006)

Model 117-04 and 177-04 grey aluminium powder coated (RAL 9007)



Basic equipment:
Head part adjustment by
gas spring

All pictures:
Model 177-04
in basic equipment;
upholstery colour 64-153
platin



Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Cut-out for echocardiography and upholstery part to close (in patient supine position left)
- Possibility to pin on a paper roll holder on both traverse sides
- Dewert-hygiene-equipment



Upholstery:

Lying surface 2-sectioned:

- Head part stepless upwards (+55°) adjustable by gas spring
- Length head part: 800 mm
- Cut-out for echocardiography and upholstery part to close, in patient supine position left, measures: L. 300 mm x W. 250 mm
- Fixed foot part
- Length foot part: 1150 mm
- Upholstery thickness: 80 mm

Technical data:

- Length of lying surface: 1950 mm
- Width of lying surface: Model 117-...: 650 mm | Model 177-...: 800 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load



Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors,
Ø 100 mm, each with total
lock




Art. no. 1140-10
Movable on 4 castors,
Ø 125 mm, each with
total lock



Art. no. 1180N
Movable on twin castors,
Ø 100 mm with central lock
(not in case of chromium-
plated frame)



Art. no. 1185
Movable on twin castors,
Ø 125 mm with central lock
(not in case of chromium-plated
frame)



Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width
650 mm / 800 mm

Upholstery options

Art. no. 1765 / 1780
Soft covering for pleasant resting
comfort, for width 650 mm / 800 mm,
thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam
cover), for width 650 mm / 800 mm,
thickness 70 mm

Side guards



Art. no. 1080-50
1 side guard lateral lowerable,
chromium-plated, completely made of steel
tube, length: 1060 mm,
length on cut-out side: 820 mm



Art. no. 1040
1 side guard tiltable, chromium-plated,
completely made of steel tube,
length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with
attachment clamp, chromium-plated,
length: 800 mm

Add-on elements



Art. no. 1410
1 piece of side rail,
length: 350 mm

Art. no. 1420
1 piece of side rail,
length: 1300 mm



Art. no. 1430
1 attachment clamp for
side rail, suitable for pipe:
Ø 18 mm



Art. no. 1440
1 arm rest for side rail,
without attachment clamp



Art. no. 1450
Infusion pole, without
attachment clamp



Art. no. 1850
Storage basket with mounting,
measures:
L. 530 x W. 370 x H. 80 mm