

Tables for echocardiography

2-sectioned table for echocardiography with hydraulic height adjustment ... for examination, treatment ...

Model 2550-00XL/H and 2551-00XL/H white powder coated (RAL 9010)
 Model 2550-03XL/H and 2551-03XL/H white aluminium powder coated (RAL 9006)
 Model 2550-04XL/H and 2551-04XL/H grey aluminium powder coated (RAL 9007)



All pictures:
 Model 2550-00XL/H
 with options: Movable on twin castors with central lock (1180N), paper roll holder (1270K);
 upholstery colour 64-708 cherry



Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- In head part: Cut-out for echocardiography and upholstery part to close (in patient supine position left)
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment

Upholstery:

Lying surface 2-sectioned:

- Head part stepless upwards (+50°) adjustable by gas spring
- Length head part: 800 mm
- In head part: Cut-out for echocardiography and upholstery part to close, in patient supine position left, measures: L. 300 mm x W. 250 mm
- Fixed foot part
- Length foot part: 1150 mm
- Upholstery thickness: 60 mm





Technical data:

- Manual foot operated height adjustment: 520 mm to 940 mm by hydraulic pump, pressure capacity of 10000 N
- Width of lying surface: Model 2550-..XL/H: 700 mm | Model 2551-..XL/H: 800 mm
- Length of lying surface: 1950 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120N-10
 Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 1180N
 Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 1185
 Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 1045N
 Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering



Paper roll holder



Art. no. 1270K / 1280K
 Paper roll holder, mounted under head part (not protruding), for width 700 mm / 800 mm



Art. no. 1270F / 1280F
 Paper roll holder, mounted at foot part (protruding: 70 mm), for width 700 mm / 800 mm

Add-on elements



Art. no. 1410
 1 piece of side rail, length: 350 mm



Art. no. 1850
 Storage basket

Side guards



Art. no. 1080-50
 1 side guard lateral lowerable, chromium-plated, completely made of steel tube, length: 1060 mm, length on cut-out side: 820 mm



Art. no. 1040
 1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
 1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

Wall protection bumpers



Art. no. 1170
 1 wall protection bumper Ø 150 mm

Upholstery options

Art. no. 1770 / 1780
 Soft covering for pleasant resting comfort, for width 700 mm / 800 mm, thickness 80 mm

Art. no. 1772 / 1782
 Upholstery medium-soft (cold foam cover), for width 700 mm / 800 mm, thickness 70 mm