

Tilt table

Tilt table with fixed height and electric motorized tilt function

... for examination, treatment, therapy ...

Model 90-00E white powder coated (RAL 9010)

Model 90-03E white aluminium powder coated (RAL 9006)

Model 90-04E grey aluminium powder coated (RAL 9007)

► Model 900-..E corresponds with model 90-..E



Model 90-00E
with options: Head part manual
upwards adjustable
(2090-058), paper roll holder
(1272), angle gauge (90-090)

Basic equipment:

- Fixed height, not adjustable
- Lying surface 1-sectioned
- Operation of electric motorized tilt function by one hand switch and one foot switch
- Integrated safety function for automatic blocking of electromotive drive
- Movable on 4 twin castors Ø 75 mm, each with lock
- Foot plate to stand on (upholstered), removable
- Side rails along lying surface (left and right side) for special tilt table equipment
- Dewert-hygiene-equipment



Upholstery:

Lying surface 1-sectioned (different head part options possible):

- Length: 1.930 mm
- Width: 700 mm
- Thickness of upholstery: 60 mm

Technical data:


- Electric motorized tilt function from -15° to +85° in 22 sec. by AC-Motor with a pressure capacity of 10000 N
- Width of lying surface: 700 mm | Width overall: 850 mm
- Length of lying surface: 1930 mm | Length overall: 1960 mm
- Load capacity: 180 kg; according to IEC 60601 tested with fourfold static load




Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Attention: necessary room height: 2.15 m



Available as an option (see pricelist for additional options)

Additional motor control options



Art. no. 1920E
Hand switch (low voltage),
with magnet on the back to
fasten



Art. no. 1910E
1 foot switch (low voltage)



Art. no. 2090-043E
Foot plate to stand on up- and
downwards adjustable by
electric motor (not removable)

Additional safety options



Art. no. 1997E
Foot switch cover,
made of brushed VA-steel

Grip handles



Art. no. 90-028
1 pair grip handles, small




Art. no. 90-038
1 pair grip handles, big

Belt elements



Art. no. 90-034
Upholstery for strap



Art. no. 90-035
1 safety belt with velcro fastener



Art. no. 90-0351
1 safety belt with automatic
lock



Art. no. 90-036
1 pair belt guides

Arm rest



Art. no. 2090-0143
Attachment clamp to fix
arm rest at side rail



Art. no. 2090-059
Arm rest

Work and rest desks



Art. no. 90-037
Work and rest desk, up- and
downwards adjustable



Art. no. 90-039
Work and rest desk, rigid

Head part options



Art. no. 90-058
Head part manual upwards
adjustable




Art. no. 90-041
Head part up- and downwards
adjustable by gas spring



Art. no. 1272
Paper roll holder, mounted under
head part (not protruding)

Special equipment



Art. no. 90-25
Extension rail



Art. no. 90-090
Angle gauge, on both sides
of the table

Upholstery options

Art. no. 90-080
Width of upholstery 800 mm