

MULTIPURPOSE TABLES  
FOR EXAMINATION, REKTOSCOPY UND GYNAECOLOGY



# SERIES 2000XL

2023.2

... WITH 3-SECTIONED LYING SURFACE




... WITH 4-SECTIONED LYING SURFACE



... WITH INCLINATION ADJUSTMENT



TABLES FOR SURGERY AND HOSPITAL



# Series 2000XL - Tables for gynaecology and rectoscopy

## Good reasons for a multipurpose table of series 2000XL

### Content

Good reasons for a multipurpose table of series 2000XL ..... 3

Multipurpose tables for examination, rectoscopy and gynaecology

4-sectioned multipurpose tables:

- Model 2640XLE ..... 6
- Model 2640XL/H ..... 8

3-sectioned multipurpose tables:

- Model 2645XLE ..... 10
- Model 2645XL/H ..... 12

3-sectioned multipurpose tables with inclination adjustment:

- Model 2648XLE ..... 14



**NEW**

All models are available in white, white aluminium and grey aluminium (frame colour).\*



MODEL 2640-03XLE



MODEL 2648-00XLE




MODEL 2645-00XLE

# Series 2000XL - Tables for gynaecology and rectoscopy

## Good reasons for a multipurpose table of series 2000XL

### Electric motorized height adjustment



- Electronical height adjustment by AC-motor with a pressure capacity of 10000 N
- Constant adjustment speed even in case of high load
- EPR-safety power cable with strain relief, length: 4 m

### Integrated safety devices (according to the newest advice of BfArM)



- Automatic blocking of electromotive drive
- Active freewheeling downwards, automatic declutch in case of resistance



- Load capacity: 220 kg



### Hydraulic height adjustment



- Manual foot operated height adjustment by a hydraulic pump with a pressure capacity of 10000 N
- Pedals for foot-controlled height positioning on both sides of the table
- Load capacity: 220 kg



### Dewert-hygiene-equipment

Enables optimal cleaning and disinfection

- Artificial leather underneath the upholstery
- Vent holes underneath the upholstery, closed by special caps with valve function
- Easy hygiene because of an open construction: good accessibility of components, no plastic coverings which are breakable in the ageing process
- Option: Innovative artificial leather Skai® Toronto EN with staynu





# Series 2000XL - Tables for gynaecology and rectoscopy

## Good reasons for a multipurpose table of series 2000XL




**Durable** ✓



**Robust** ✓




**Easy-care** ✓



**Twin frame  
scissor lifting system** ✓



**Proven  
Dewert-quality**




**High  
adjustment speed** ✓

### Technology:



- Extremely stable twin frame scissor lifting system, with high torsional stiffness due to thick-walled square tubes with big dimensions, also for under and upper frame and scissor

- High stability due to massive 20 mm axes



- Widely-spaced feet and a high deadweight guarantee stability

- Robust frame construction which is proven in demanding everyday life of hospitals for years



- Abrasion-resistant guide system at under and upper frame due to the use of roller bearings



### Upholstery:

- Thickness of upholstery: 60 mm
- Multiple glued, high-strength laminated wood panel (thickness: 18 mm) no breaking chipboards in case of pressure on the edges

- Two-layer homogenous polyurethane foam, density 46 and 75, Ökotex-certified

- 2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Cleaning kit for Skai products  
(Art. no. 050-1000)

Consists of:

- Spray bottle with 200 ml special Skai cleansing agent
- 1 microfibre care cloth 40 x 40 cm

# Series 2000XL - Tables for gynaecology and rectoscopy

## Good reasons for a multipurpose table of series 2000XL

### ... with 3-sectioned or 4-sectioned lying surface:

- Head part 30° positive adjustable by 2 metal ratchets
- Middle part 50° positive adjustable by rack system
- Seat section with half-round cut-out, manually inclinable to +20°
- 4-sectioned lying surface additionally: Rectoscopy insert to lift up, width: 100 mm



- 3-sectioned lying surface and inclination adjustment additionally: Electric motorized inclination adjustment to -15°, operation by hand switch



### Equipment:



- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat basin, support: HPL panel
- Movable on 4 comfort twin castors, Ø 100 mm, each with lock
- Clamps for leg supports, variants of leg supports available
- Dewert-hygiene-equipment



You will find more options on the following pages.




# Series 2000XL

## 4-sectioned multipurpose table with electric motorized height adjustment ... for examination, rectoscopy, gynaecology ...

**Model 2640-00XLE white powder coated (RAL 9010)**

**Model 2640-03XLE white aluminium powder coated (RAL 9006)**

**Model 2640-04XLE grey aluminium powder coated (RAL 9007)**



All pictures::

Model 2640-03XLE with options:  
loss of basic mobility (1182),  
movable on twin castors with central lock (1180N),  
head part adjustable by gas spring (1655), head part with  
nose opening and upholstery part to close (1610, 1620),  
paper roll holder (1265K), side rail (1410)  
upper picture additionally: Leg supports acc. to  
Goepel (1520); Upholstery colour: 64-701 stone blue

### Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Electromotive drive control by foot switch (option: hand switch)
- Integrated safety function for automatic blocking of electromotive drive (art. no. 1999, mandatory equipment)
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment



### Upholstery:

Lying surface 4-sectioned:

- Head part 30° positive adjustable by 2 metal rachets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Rectoscopy insert to lift up and fixable, width: 100 mm
- Seat section with half-round cut-out, manually inclinable (20°)
- Length seat section: 400 mm
- Thickness of upholstery: 80 mm



### Technical data:

- Electronical height adjustment: 620 mm to 1060 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1980 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



### Design:


2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



**Available as an option (see pricelist for additional options)**

**Mobility**



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm

**Nose opening**



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close



**Additional motor control options**



Art. no. 1920E  
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 1910E  
Foot switch (low voltage)



Art. no. 1940E  
1 foot switch fixation mounted at underframe



Art. no. 1950E  
Foot switch rail, lateral, on both sides, security adjustment

**Paper roll holder**



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for with of lying surface 700 mm

**Side guards**



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

**Head part options**



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

**Add-on elements**



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

**Leg supports**



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

**Upholstery options**

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

# Series 2000XL


## 4-sectioned multipurpose table with hydraulic height adjustment

... for examination, rectoscopy, gynaecology ...

**Model 2640-00XL/H white powder coated (RAL 9010)**

**Model 2640-03XL/H white aluminium powder coated (RAL 9006)**

**Model 2640-04XL/H grey aluminium powder coated (RAL 9007)**



2 upper pictures:  
Model 2640-00XL/H  
with option:  
Paper roll holder (1265K),  
side rail (1410),  
leg supports (1510);  
upholstery colour: 64-174 sky

(Model -01, top frame  
chromium-plated  
no longer available)



### Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment



### Upholstery:

Lying surface 4-sectioned:

- Head part 30° positive adjustable by 2 metal ratchets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Rectoscopy insert to lift up and fixable, width: 100 mm
- Seat section with half-round cut-out, manually inclinable (20°)
- Length seat section: 400 mm
- Thickness of upholstery: 80 mm



### Technical data:

- Manual foot operated height adjustment: 640 mm to 1060 mm by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1980 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

### Mobility



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm

### Nose opening



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close

### Paper roll holder



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for width of lying surface 700 mm

### Side guards



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

### Head part options



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

### Add-on elements



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

### Leg supports



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

### Upholstery options

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

# Series 2000XL


## 3-sectioned multipurpose table with electric motorized height adjustment

... for examination, rectoscopy, gynaecology ...

**Model 2645-00XLE white powder coated (RAL 9010)**

**Model 2645-03XLE white aluminium powder coated (RAL 9006)**

**Model 2645-04XLE grey aluminium powder coated (RAL 9007)**



Model 2645-00XLE  
with option :  
Movable on twin castors  
with central lock (1180N);  
upper picture additionally:  
Leg supports (1510);  
Upholstery colour:  
64-173 smaragd



### Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Electromotive drive control by foot switch (option: hand switch)
- Integrated safety function for automatic blocking of electromotive drive (art. no. 1999, mandatory equipment)
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment



### Upholstery:

Lying surface 3-teilig:

- Head part 30° positive adjustable by 2 metal rachets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Seat section with half-round cut-out, manually inclinable (20°)
- Length seat section: 450 mm
- Thickness of upholstery: 80 mm



### Technical data:

- Electronical height adjustment: 620 mm to 1060 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1980 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



### Design:


2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

**Mobility**



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm

**Nose opening**



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close



**Additional motor control options**



Art. no. 1920E  
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 1910E  
Foot switch (low voltage)



Art. no. 1940E  
1 foot switch fixation mounted at underframe



Art. no. 1950E  
Foot switch rail, lateral, on both sides, security adjustment

**Paper roll holder**



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for with of lying surface 700 mm

**Side guards**



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

**Head part options**



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

**Add-on elements**



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

**Leg supports**



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

**Upholstery options**

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm



# Series 2000XL


## 3-sectioned multipurpose table with hydraulic height adjustment

... for examination, rectoscopy, gynaecology ...

Model 2645-00XL/H white powder coated (RAL 9010)

Model 2645-03XL/H white aluminium powder coated (RAL 9006)

Model 2645-04XL/H grey aluminium powder coated (RAL 9007)



All pictures:  
Model 2645-01XL/H  
with option:  
Paper roll holder (1265K);  
upper 2 pictures additionally:  
Leg supports (1510);  
upholstery colour: 64-717 deep sea

(Model -01, top frame  
chromium-plated no longer  
available)

### Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment



### Upholstery:

Lying surface 3-sectioned:

- Head part 30° positive adjustable by 2 metal ratchets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Seat section with half-round cutout, manually inclinable (20°)
- Length seat section: 450 mm
- Thickness of upholstery: 80 mm

### Technical data:

- Manual foot operated height adjustment: 640 mm to 1060 mm by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1980 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)


### Mobility



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



### Nose opening



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close

### Paper roll holder



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for width of lying surface 700 mm

### Side guards



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

### Head part options



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

### Add-on elements



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

### Leg supports



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

### Upholstery options

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

# Series 2000XL


## 3-sectioned multipurpose table with electric motorized height and inclination adjustment

... for examination, rectoscopy, gynaecology ...

**Model 2648-00XLE white powder coated (RAL 9010)**

**Model 2648-03XLE white aluminium powder coated (RAL 9006)**

**Model 2648-04XLE grey aluminium powder coated (RAL 9007)**



Model 2648-00XLE with options:  
Loss of basic mobility (1182); movable on twin castors with central lock (1180N), paper roll holder (1265K); upper picture additionally:  
Leg supports acc. to Goepel (1520); Upholstery colour: 64-715 anthracite



### Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Electromotive drive control of height adjustment by foot switch
- Electromotive drive control of inclination adjustment by hand switch, inclination adjustment to -15°
- Integrated safety function for automatic blocking of electromotive drive (art. no. 1999, mandatory equipment)
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment



### Upholstery:

Lying surface 3-sectioned:

- Head part 30° positive adjustable by 2 metal rachets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Seat section with half-round cutout, manually inclinable (20°)
- Length seat section: 450 mm
- Thickness of upholstery: 80 mm



### Technical data:

- Electronical height adjustment: 650 mm to 1090 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1950 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

### Mobility



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm

### Nose opening



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close



### Additional motor control options



Art. no. 1920E  
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 1910E  
Foot switch (low voltage)



Art. no. 1940E  
1 foot switch fixation mounted at underframe



Art. no. 1950E  
Foot switch rail, lateral, on both sides, security adjustment

### Paper roll holder



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for with of lying surface 700 mm

### Side guards



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

### Head part options



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

### Add-on elements



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

### Leg supports



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

### Upholstery options

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm



Technical data subject to changes and subject to errors. Due to technical reasons colour deviations of shown and origin colours are possible. Products are not approved for the American market.

K.H. Dewert GmbH  
Vollmestraße 7  
33649 Bielefeld  
Germany

Phone + 49 (0) 5 21 / 4 00 27-0  
Fax + 49 (0) 5 21 / 4 00 27-27

Internet: [www.khdewert.de](http://www.khdewert.de)  
E-Mail: [info@khdewert.de](mailto:info@khdewert.de)  
YouTube: K.H. Dewert

# MULTIPRPOSE TABLES FOR EXAMINATION, RECTOSCOPY, GYNAECOLOGY

You will find the whole Dewert product range with additional articles and equipment on our website and in our further catalogues. We look forward to hear from you!



04.2023\_03

TABLES FOR SURGERY AND HOSPITAL | SERIES 2000XL

