Series DFH

Short tables for examination and treatment

1-sectioned table with fixed height

Model 16-00 and 17-00 white powder coated (RAL 9010)
Model 16-03 and 17-03 white aluminium powder coated (RAL 9006)
Model 16-04 and 17-04 grey aluminium powder coated (RAL 9007)

German manufacture



Both pictures:
Model 16-00 with option:
Paper roll holder (1210);
upholstery colour: 64-708 cherry

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Possibility to pin on a paper roll holder on both traverse sides
- Dewert-hygiene-equipment



Upholstery:

Lying surface 1-sectioned:

- Length: 1450 mm (other measures possible on request)
- Upholstery thickness: 60 mm

Technical data:

- Length overall: 1450 mm
- Width overall: Model 16-..: 650 mm | Model 17-..: 800 mm
- Height: 650 mm
- Load capacity: 250 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-10
Movable on 4 castors, Ø 125 mm, each with total lock



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock



Art. no. 1185
Movable on twin castors, Ø 125 mm with central lock

Longitudinal brace



Art. no. 1310
Longitudinal brace for additional stability, necessary in case of mobility

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam cover), for width 650 mm / 800 mm, thickness 70 mm

DEWERT Präzisions-Stahlrohrmöbel

Updated 06.2022_01
changes concerning the technical data are subject to change