Series DFH

2-sectioned table with fixed height made of square tube
Length head part: 750 mm

Model 20-00 white powder coated (RAL 9010)
Model 20-03 white aluminium powder coated (RAL 9006)
Model 20-04 grey aluminium powder coated (RAL 9007)

German manufacture



Both pictures:
Model 20-04
with option:
Paper roll holder (1210);
upholstery colour: 64-708 cherry

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Head part stepless adjustable to sitting position
- Possibility to pin on a paper roll holder on head side and foot side
- Dewert-hygiene-equipment

Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable to +75° by gas spring (sitting position)
- Length head part: 750 mm
- Fixed foot part
- Length foot part: 1200 mm
- Upholstery thickness: 60 mm



Technical data:

- Length overall: 1950 mm
- Width overall: 650 mm
- Height: 680 mm
- Load capacity: 250 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-10
Movable on 4 castors, Ø 125 mm, each with total lock



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock



Art. no. 1185
Movable on twin castors, Ø 125 mm with central lock

Longitudinal brace



Art. no. 1310
Longitudinal brace for additional stability, necessary in case of mobility

Head part options



Art. no. 1662
Additional second head part of the same size, adjustable by 2 metal rachets



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Side guards



Art. no. 1088
2 side guards lateral lowerable, RAL 9006 powder coated, length: 820 mm



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm,

Paper roll holder



Art. no. 1210
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765
Soft covering for pleasant resting comfort, for width 650 mm, thickness 80 mm

Art. no. 1762
Upholstery medium-soft (cold foam cover), for width 650 mm, thickness 70 mm



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm

Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm

Art. no. 1420
1 piece of side rail, length: 1300 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1440
1 arm rest for side rail, without attachment clamp



Art. no. 1450
Infusion pole, without attachment clamp

Updated 06.2022_01
changes concerning the technical data are subject to change